



CooperVision®
Live Brightly.®

Product Specifications Guide



MyDay® daily disposable family

MyDay® daily disposable



High-performing 1 Day contact lenses for demanding lifestyles.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
stenfilcon A / 54%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

Sphere Power
-0.25D to -6.00D (0.25D steps)
-6.50D to -12.00D (0.50D steps)
*+0.25D to +5.00D (0.25D steps)
*+5.50D to +6.00D (0.50D steps)

Design
Aspheric

MyDay® Toric



High-performing 1 Day contact lenses for demanding lifestyles.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
Stenfilcon A / 54%

Base Curve 8.6mm
Diameter 14.5mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes
Design
Aspheric

Sphere Power
Plano to -6.00D (0.25D steps)
-6.50D to -10.00D (0.50D steps)
*+0.50D to +6.00D (0.25D steps)

Cylinder Power
-0.75
-1.25
-1.75
-2.25

Axis
10° to 180° (10° steps)
10°, 20°, 90°, 160°, 170°, 180°

-0.75
-1.25
-1.75
-2.25
10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°
10°, 20°, 90°, 160°, 170°, 180°

MyDay® Multifocal



High-performing 1 Day contact lenses for demanding lifestyles.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
stenfilcon A / 54%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

Sphere Power
+4.00D to -10.00D (0.25D steps)
-10.50D to -12.00D (0.50D steps)
*+4.25D to +8.00D (0.25D steps)

Addition:
Low Addition +0.75D to +1.25D
Medium Addition +1.50D to +1.75D
High Addition +2.00D to +2.50D



CooperVision
Live Brightly!

To learn more about our contact lenses visit
cooper-vision.com/contact-lenses

clariti® 1 day family

clariti® 1 day



High breathability¹, comfort and great value[†] in convenient daily disposable contact lenses.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
somofilcon A / 56%

Base Curve 8.6mm
Diameter 14.1mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

Sphere Power
no plus power:
-0.50D to -6.00D
(0.25D steps)

Design
Aspheric

-6.50D to -10.00D
(0.50D steps)

clariti® 1 day toric



High breathability¹, comfort and great value[†] in convenient daily disposable contact lenses.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
somofilcon A / 56%

Base Curve 8.6mm
Diameter 14.3mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

Sphere Power
Plano to -6.00D
(0.25D steps)

Cylinder Power
-0.75
-1.25
-1.75
-2.25

Axis
10° to 180°
(10 steps)

10°, 20°, 70°, 80°,
90°, 100°, 110°,
160°, 170°, 180°

Design
Aspheric

-6.50D to -9.00D
(0.50D steps)

-0.75
-1.25
-1.75
-2.25

10°, 20°, 60°, 70°,
80°, 90°, 100°, 110°,
120°, 160°, 170°, 180°

10°, 20°, 90°,
160°, 170°, 180°

*+0.25D to +4.00D
(0.25D steps)

-0.75
-1.25
-1.75

10°, 20°, 70°, 80°,
90°, 100°, 110°,
160°, 170°, 180°

clariti® 1 day multifocal



High breathability¹, comfort and great value[†] in convenient daily disposable contact lenses.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
somofilcon A / 56%

Base Curve 8.6mm
Diameter 14.1mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint No

Sphere Power
+5.00D to -6.00D
(0.25D steps)

Add Power
Low:
up to +2.25D

High:
+2.50D to
+3.00D

Design
center near,
progressive
intermediate
and peripheral
distance

Biofinity® family

Biofinity®



Wearing Schedule

Daily or up to 29 nights continuous

Carton Size Options

6 pack

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.6mm

Diameter 14.0mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Sphere Power

+8.00D to -12.00D
(0.50D steps after +/-6.00D)

Design

Aspheric

Premium monthly replacement contact lenses for long-lasting comfort and excellent vision²⁻⁴

Biofinity® toric



Wearing Schedule

Daily or up to 29 nights continuous

Carton Size Options

6 pack

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.7mm

Diameter 14.5mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Sphere Power

+8.00D to -10.00D
(0.50D steps after +/-6.00D)

Cylinder Power

-0.75
-1.25
-1.75
-2.25

Axis

10° to 180°
(10° steps)

Design

Optimised Toric Lens Geometry™

Premium monthly replacement contact lenses for long-lasting comfort and excellent vision²⁻⁴

Biofinity® multifocal



Wearing Schedule

Daily or up to 29 nights continuous

Carton Size Options

6 pack

Material / H₂O Content

comfilcon A / 48%

Base Curve 8.6mm

Diameter 14.0mm

FDA Group SiHy

UV Blocking* No

Visibility Tint Yes

Sphere Power

+6.00D to -10.00D
(0.50D steps after -6.00D)

Add Power

+1.00
+1.50
+2.00
+2.50

Design

D lens
N lens

Premium monthly replacement contact lenses for presbyopia, to provide you incredible comfort¹ and excellent vision at all distances⁵




CooperVision[®]
Live Brightly![™]

To learn more about our contact lenses visit
cooper-vision.com/contact-lenses

Avaira Vitality[™] family

Avaira Vitality[™]



Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
fanfilcon A / 55%

Base Curve 8.4mm
Diameter 14.2mm
FDA Group SiHy
UV Blocking* Yes
Visibility Tint Yes

Sphere Power
-0.25D to -6.00D
(0.25D steps)

-6.50D to -12.00D
(0.50D steps)

Design
Aspheric front surface

Contact lenses with UV blocking:
everything you expect from
reusable contact lenses.

MiSight[®] 1 day Myopia Management

MiSight[®] 1 day Myopia Management



Wearing Schedule
Daily

Carton Size Options
30 pack

Material / H₂O Content
omafilcon A

Base Curve 8.7mm
Diameter 14.2mm
FDA Group 2
UV Blocking* No
Visibility Tint Yes

Sphere Power
-0.25D to -6.00D
(0.25D steps)

-6.50D to -10.00D
(0.50D steps)

Child-friendly daily disposable
contact lenses proven to slow the
progression of nearsightedness in
children.⁶



CooperVision®
Live Brightly!

To learn more about our contact lenses visit
cooper-vision.com/contact-lenses

Proclear® family

Proclear® 1 day



Clear vision*, convenience and comfort, for contact lens wearers who experience eye dryness.

Wearing Schedule
Daily

CARTON SIZE OPTIONS
30 pack

Material / H₂O Content
omafilcon A / 60%

Base Curve 8.7mm
Diameter 14.2mm
FDA Group 2

Sphere Power
-0.25D to -6.00D
(0.25D steps)

Design
Aspheric

UV Blocking* No
Visibility Tint Yes

-6.50D to -12.00D
(0.50D steps)

*+0.50D to +6.00D
(0.25D steps)

+6.50D to +8.00D
(0.50D steps)

Proclear® 1 day multifocal



Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

Wearing Schedule
Daily disposable

Carton Size Options
30 pack

Material / H₂O Content
omafilcon A / 60%

Base Curve 8.7mm
Diameter 14.2mm
FDA Group 2

Sphere Power
Plano to -10.00D
(0.50D steps after -6.00D)

Add Power
Single power profile

Design
Center near
Aspheric Design

UV Blocking* No
Visibility Tint Yes

+0.25D to +6.00D
(0.25D steps)

Proclear® multifocal



Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

Wearing Schedule
Daily

Carton Size Options
6 pack

Material / H₂O Content
omafilcon B / 62%

Base Curve 8.7mm
Diameter 14.4mm
FDA Group 2

Sphere Power
-0.50D to -10.00D
(0.50D steps after -6.00D)

Add Power
+1.00
+1.50
+2.00
+2.50

Design
D lens
N lens

UV Blocking* No
Visibility Tint Yes



CooperVision®
Live Brightly!

To learn more about our contact lenses visit
cooper-vision.com/contact-lenses

Biomedics® family

Biomedics® 1 day Extra



Wearing Schedule
Daily

Carton Size Options
30 pack

Material / H₂O Content
ocufilcon D / 55%

Base Curve 8.6mm
Diameter 14.2mm

Sphere Power
-0.25D to -10.00D
(0.50D steps after -6.00D)

Design
Aspheric

Designed to maximise the feeling of natural all-day comfort.

Biomedics® 1 day Extra toric



Wearing Schedule
Daily

Carton Size Options
30 pack

Material / H₂O Content
ocufilcon D / 55%

Base Curve 8.7mm
Diameter 14.5mm

Sphere Power
Plano to -6.00D
(0.25D steps)

Cylinder Power
-0.75
-1.25
-1.75

Axis
20°, 90°,
160°, 180°

All the great things about daily disposable for those with astigmatism.

-6.50D to -10.00D
(0.50D steps)

-0.75
-1.25
-1.75,

90°, 180°

¹ Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005. 2005;82(6):467-472.

² CooperVision data on file. 2006. Aquaform® Technology uses a matrix of long silicone chains with hydrogen bonds to lock water molecules in the lens for incredible, long-lasting comfort and clarity.

³ CooperVision data on file. 2014. Clinical evaluation of Biofinity® spherical contact lens compared to Air Optix® Aqua spherical lenses over one month of wear.

⁴ CooperVision data on file. 2014. Clinical evaluation of Biofinity® spherical contact lens compared to Ultra spherical contact lens over one month of wear.

⁵ CooperVision data on file 2011. Clinical Evaluation of Biofinity multifocal® compared to ACUVUE® Oasys for Presbyopia, AIR OPTIX® for Presbyopia and PureVision® Multi-Focal.

⁶ Chamberlain P, et al. A 3-year randomized clinical trial of MiSight® lenses for myopia control. Optom Vis Sci. 2019;96(8):556-567.

⁷ CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

*UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.

All trademarks, service marks, trade names, trade dress and product names featured herein are registered trademarks and trademarks of CooperVision, its affiliates and subsidiaries, unless otherwise noted. CooperVision's trademarks may be used publicly only with permission.