



CooperVision®  
Live Brightly.®

# Product Specifications Guide



## MyDay® daily disposable family

### MyDay® daily disposable



High-performing 1 Day contact lenses for demanding lifestyles.

**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
stenfilcon A / 54%

**Base Curve** 8.4mm  
**Diameter** 14.2mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** Yes

**Sphere Power**  
-0.50 to -6.00D  
(0.25D steps)

**Design**  
Aspheric

-6.50 to -10.00D  
(0.50D steps)

## clariti® 1 day family

### clariti® 1 day



High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
somofilcon A / 56%

**Base Curve** 8.6mm  
**Diameter** 14.1mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** No

**Sphere Power**  
-0.50D to -10.00D  
+0.50D to +8.00D  
(0.50D steps after  
+/-6.00D)  
No Plano

**Design**  
Aspheric

### clariti® 1 day toric



High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
somofilcon A / 56%

**Base Curve** 8.6mm  
**Diameter** 14.3mm  
**FDA Group** SiHy  
**UV Blocking\*** Yes  
**Visibility Tint** No  
**Design**  
Aspheric

**Sphere Power**  
Plano to - 9.00D  
(0.50D steps after  
-6.00D)

**Cylinder Power**  
-0.75  
-1.25  
-1.75

**Axis**  
10°, 20°, 60°,  
70°, 80°, 90°,  
100°, 110°, 120°,  
160°, 170°, 180°

+0.25D to +4.00D  
(0.25D steps)

-0.75  
-1.25  
-1.75

20°, 70°, 90°,  
110°, 160°, 180°

Plano to -9.00D  
(0.50D steps after  
-6.00D)

-2.25

0°, 20°, 90°,  
160°, 170°, 180°



CooperVision<sup>®</sup>  
Live Brightly.<sup>™</sup>

To learn more about our contact lenses visit  
[coopervision.com/contact-lenses](http://coopervision.com/contact-lenses)

## clariti<sup>®</sup> 1 day multifocal



High breathability<sup>1</sup>, comfort and great value<sup>†</sup> in convenient daily disposable contact lenses.

### Wearing Schedule

Daily disposable

### Carton Size Options

30 pack

### Material / H<sub>2</sub>O Content

somofilcon A / 56%

**Base Curve** 8.6mm

**Diameter** 14.1mm

**FDA Group** SiHy

**UV Blocking\*** Yes

**Visibility Tint** No

### Sphere Power

+5.00D to -6.00D  
(0.25D steps)

### Add Power

Low:  
up to +2.25D

High:  
+2.50D to  
+3.00D

### Design

center near,  
progressive  
intermediate  
and peripheral  
distance

## Biofinity<sup>®</sup> family

### Biofinity<sup>®</sup>



Premium monthly replacement contact lenses for long-lasting comfort and excellent vision<sup>2-4</sup>

### Wearing Schedule

Daily or up to 29  
nights continuous

### Carton Size Options

6 pack

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

### Sphere Power

+8.00D to -12.00D  
(0.50 steps after  
+/-6.00D)

### Design

Aspheric

### Biofinity<sup>®</sup> toric



Premium monthly replacement contact lenses for long-lasting comfort and excellent vision<sup>2-4</sup>

### Wearing Schedule

Daily or up to 29  
nights continuous

### Carton Size Options

6 pack

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

**Base Curve** 8.7mm

**Diameter** 14.5mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

### Sphere Power

+8.00D to -10.00D  
(0.50 steps after  
+/-6.00D)

### Cylinder Power

-0.75  
-1.25  
-1.75  
-2.25

### Axis

10° to 180°  
(10° steps)

### Design

Optimised Toric  
Lens Geometry<sup>™</sup>



CooperVision<sup>®</sup>  
Live Brightly!<sup>™</sup>

To learn more about our contact lenses visit  
[coopervision.com/contact-lenses](http://coopervision.com/contact-lenses)

## Biofinity<sup>®</sup> multifocal



### Wearing Schedule

Daily or up to 29 nights continuous

### Sphere Power

+6.00D to -10.00D  
(0.50D steps after -6.00D)

### Carton Size Options

3 pack

### Add Power

+1.00  
+1.50  
+2.00  
+2.50

### Material / H<sub>2</sub>O Content

comfilcon A / 48%

### Design

D lens  
N lens

**Base Curve** 8.6mm

**Diameter** 14.0mm

**FDA Group** SiHy

**UV Blocking\*** No

**Visibility Tint** Yes

Premium monthly replacement contact lenses for presbyopia, to provide you incredible comfort<sup>1</sup> and excellent vision at all distances<sup>5</sup>

## Avaira Vitality<sup>™</sup> family

### Avaira Vitality<sup>™</sup>



### Wearing Schedule

Daily

### Sphere Power

-0.50 to -6.00DS  
(0.25D steps)

-6.50 to -12.00DS  
(0.50D steps)

### Carton Size Options

6 pack

### Design

Aspheric front surface

### Material / H<sub>2</sub>O Content

fanfilcon A / 55%

**Base Curve** 8.4mm

**Diameter** 14.2mm

**FDA Group** SiHy


**UV Blocking\*** Yes

**Visibility Tint** Yes

Contact lenses with UV blocking: everything you expect from reusable contact lenses.

## Brilliant Futures<sup>™</sup> with MiSight<sup>®</sup> Myopia Management

### Brilliant Futures<sup>™</sup> with MiSight<sup>®</sup> Myopia Management



### Wearing Schedule

Daily

### Sphere Power

-0.25 to -6.00  
(0.25 steps)

### Carton Size Options

30 pack

### Material / H<sub>2</sub>O Content

omafilcon A

**Base Curve** 8.7mm

**Diameter** 14.2mm

**FDA Group** 2

**UV Blocking\*** No

**Visibility Tint** Yes

Child-friendly daily disposable contact lenses proven to slow the progression of nearsightedness in children.<sup>6</sup>



CooperVision<sup>®</sup>  
Live Brightly!<sup>™</sup>

To learn more about our contact lenses visit  
[coopervision.com/contact-lenses](https://coopervision.com/contact-lenses)

## Proclear<sup>®</sup> family

### Proclear<sup>®</sup> 1 day



Clear vision\*, convenience and comfort, for contact lens wearers who experience eye dryness.

**Wearing Schedule**  
Daily

**CARTON SIZE OPTIONS**  
30 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A / 60%

**Base Curve** 8.7mm  
**Diameter** 14.2mm  
**FDA Group** 2  
**UV Blocking\*** No  
**Visibility Tint** Yes

**Sphere Power**  
-0.50 to -10.00D  
(0.50 steps after -6.00D)

**Design**  
Aspheric

+0.25 to +8.00  
(0.50 steps after +6.00D)

### Proclear<sup>®</sup> 1 day multifocal



Convenience and comfort for multifocal contact lens wearers who experience eye dryness.

**Wearing Schedule**  
Daily disposable

**Carton Size Options**  
30 pack

**Material / H<sub>2</sub>O Content**  
omafilcon A / 60%

**Base Curve** 8.7mm  
**Diameter** 14.2mm  
**FDA Group** 2  
**UV Blocking\*** No  
**Visibility Tint** Yes

**Sphere Power**  
Plano to -10.00D  
(0.50 steps after -6.00)

**Add Power**  
Single power profile

**Design**  
Center near  
Aspheric Design

+0.25 to +6.00  
(0.25 steps)

### Proclear<sup>®</sup>



Clear vision<sup>7</sup> and comfort for contact lens wearers who experience eye dryness.

**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.6mm  
**Diameter** 14.2mm  
**FDA Group** 2  
**UV Blocking\*** No  
**Visibility Tint** Yes

**Sphere Power**  
-0.50D to -10.00D  
(0.50 steps after -6.00)

**Design**  
Sphere

## Proclear® multifocal



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
omafilcon B / 62%

**Base Curve** 8.7mm  
**Diameter** 14.4mm  
**FDA Group** 2

**Sphere Power**  
-0.50D to -10.00D  
(0.50D steps after -6.00D)

**Add Power**  
+1.00  
+1.50  
+2.00  
+2.50

**Design**  
D lens  
N lens

**UV Blocking\*** No  
**Visibility Tint** Yes

Helps provide all-day comfort for multifocal contact lens wearers who experience eye dryness.

## Biomedics® family

### Biomedics® 55 Evolution



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
ocufilcon D / 55%

**Base Curve** 8.8  
**Diameter** 14.2mm  
**FDA Group** 4

**Sphere Power**  
-0.50 to -10.00D  
(0.50 steps after -6.00)

**Design**  
Aspheric

**UV Blocking\*** No  
**Visibility Tint** Yes

Designed for comfort and performance at an affordable price.

### Biomedics® toric



**Wearing Schedule**  
Daily

**Carton Size Options**  
6 pack

**Material / H<sub>2</sub>O Content**  
ocufilcon D / 55%

**Base Curve** 8.7mm  
**Diameter** 14.5mm  
**FDA Group** 4

**Sphere Power**  
Plano to -9.00D  
(0.50 steps after -6.00D)

**Cylinder Power**  
-0.75  
-1.25  
-1.75  
-2.25

**Axis**  
10° to 180°  
(10° steps)

**UV Blocking\*** Yes  
**Visibility Tint** Yes

The affordable contact lens for wearers with astigmatism.

## Biomedics® 1 day Extra



<b>Wearing Schedule</b>	Daily	<b>Carton Size Options</b>	30 pack	<b>Material / H<sub>2</sub>O Content</b>	ocufilcon D / 55%	<b>Base Curve</b>	8.6mm
						<b>Diameter</b>	14.2mm

<b>Sphere Power</b>	-0.50D to -10.00D (0.50D steps after -6.00D)	<b>Design</b>	Aspheric
---------------------	---	---------------	----------

Designed to maximise the feeling of natural all-day comfort.

## Biomedics® 1 day Extra toric



<b>Wearing Schedule</b>	Daily	<b>Carton Size Options</b>	30 pack	<b>Material / H<sub>2</sub>O Content</b>	ocufilcon D / 55%	<b>Base Curve</b>	8.7mm
						<b>Diameter</b>	14.5mm

<b>Sphere Power</b>	-6.00 to 0.00 (0.25D steps)	<b>Cylinder Power</b>	-0.75, -1.25, -1.75	<b>Axis</b>	20°, 90°, 160°, 180°
	-7.50 to 6.00 (0.50D steps)		-0.75, -1.25, -1.75		20°, 90°, 160°, 180°
	-10.00 to -7.50 (0.50D steps)		-0.75, -1.25, -1.75		90°, 180°

All the great things about daily disposable for those with astigmatism.

<sup>1</sup> Brennan NA. Beyond Flux: Total Corneal Oxygen Consumption as an Index of Corneal Oxygen During Contact Lens Wear. Optom Vis Sci 2005;82(6):467-472.

<sup>2</sup> CooperVision data on file. 2006. Aquaform® Technology uses a matrix of long silicone chains with hydrogen bonds to lock water molecules in the lens for incredible, long-lasting comfort and clarity.

<sup>3</sup> CooperVision data on file. 2014. Clinical evaluation of Biofinity® spherical contact lens compared to Air Optix® Aqua spherical lenses over one month of wear.

<sup>4</sup> CooperVision data on file. 2014. Clinical evaluation of Biofinity® spherical contact lens compared to Ultra spherical contact lens over one month of wear.

<sup>5</sup> CooperVision data on file 2011. Clinical Evaluation of Biofinity multifocal® compared to ACUVUE® Oasys for Presbyopia, AIR OPTIX® for Presbyopia and PureVision® Multi-Focal.

<sup>6</sup> Chamberlain P, et al. A 3-year randomized clinical trial of MiSight® lenses for myopia control. Optom Vis Sci. 2019;96(8):556-567.

<sup>7</sup> CVI data on file 2017 & 2014; non-dispensing short-term fitting study. N=40 & n=20.

\*UV-blocking contact lenses help provide protection against transmission of harmful UV radiation to eye but are not substitutes for protective UV-absorbing eyewear, such as UV absorbing goggles or sunglasses, as they do not completely cover the eye or surrounding area. Continue to use UV-absorbing eyewear as directed by your eye care professional.

All trademarks, service marks, trade names, trade dress and product names featured herein are registered trademarks and trademarks of CooperVision, its affiliates and subsidiaries, unless otherwise noted. CooperVision's trademarks may be used publicly only with permission.